

Microbiology Collection / Transport

Routine aerobic bacterial culture use either Amies or Liquid Stuart transport: Amies Transport Media (Blue) Aerobic bacterial culture (MC, TC, GENC, GBBSC, GBSP, etc..)



Copan Venturi Transystem Double Swab PCR testing. (CROSP, MRSAP, VREP & SSTIP)



eSwab

Anaerobic bacterial culture (AC) – Also acceptable for aerobic bacterial culture. (MC, TC, GENC, GBBSC, etc..)



Aptima Urine Collection Kit

Transportation of amplified GC & Chlamydia on urine specimens. (CNGAM, STD3 & MGEN)



Aptima Unisex Collection Kit

Amplified GC & Chlamydia (CNGAM, STD3 & MGEN)



Aptima Multitest Swab Collection Kit

Amplified GC & Chlamydia (CNGAM, STD3 & MGEN)



Universal Transport Media (UTM) / Copan Flocked Swab - Viral collection & transport (HSVAM, VZVAM, HVZAM, COVP, RP4, INFAB, INFP, INFRP, RPPCR & RSVR)

BD Affirm ATTS

Vaginal pathogens by DNA Probes (VAGP)



Urine Culture Transport Kit

Use for transport of urine samples for bacterial culture (UC & USR)



C&S Medium (Liquid Cary Blair)

Use for stool culture and routine O&P testing. (FCSHG, FSHGC, FC, CAMAG, YER, OVP, GAT, CRYPT, GICPP & GIVPP)



UA Preservative Tube

Urine dipstick & microscopic (US, USM, UM, UMCRO, USR)



O&P Transport Kit

Use if history of travel outside the US or immunocompromised (OPFA)



Key: () - indicates test code(s) used with appropriate transport